


Vascular Intervention // Coronary
Coronary Guide Wire

Galeo Pro



Control: Enhanced torque response
for steerability



Feedback: Minimized friction
for tactile distal feedback



Choice: Broad selection of tip
flexibilities and support levels

Galeo Pro

Vascular
Intervention
Coronary



Indicated to facilitate the placement of interventional cardiology catheters.*

Technical Data		Guide wire					
	Diameter	0.014"					
	Length	190 cm and 300 cm					
	Core wire material	High tensile stainless steel					
	Proximal coil	21.5 cm, stainless steel					
	Distal coil	26 mm, palladium, radiopaque					
	Proximal (shaft) coating	PTFE					
	Distal coating	Hydrophilic (30 cm)					
	Range of tip flexibilities	High Flexible (HF), Flexible (F), Medium (M)					
	Range of support levels	Standard and Extra Support (ES)					
	Shaft markers	92 cm and 102 cm from distal end					
Ordering Information		Catalogue number	Tip shape	Tip flexibility	Label color code	Support level	Length
	389781	Straight	High flexible	Red	Standard	190 cm	
	389782	J	High flexible	Red	Standard	190 cm	
	389783	Straight	Flexible	Orange	Standard	190 cm	
	389784	J	Flexible	Orange	Standard	190 cm	
	389785	Straight	Medium	Blue	Standard	190 cm	
	389786	J	Medium	Blue	Standard	190 cm	
	389787	Straight	Flexible	Yellow	Extra support	190 cm	
	389788	J	Flexible	Yellow	Extra support	190 cm	
	389791	Straight	Flexible	Orange	Standard	300 cm	

*Indication as per IFU.

Manufacturer:
Concert Medical, LLC
77 Accord Park Dr.
Norwell, MA 02061-1605, USA
Tel +1 (781) 871 7882
info@concertmedical.com
www.concertmedical.com

Distributor:
BIOTRONIK AG
Ackerstrasse 6
8180 Bülach, Switzerland
Tel +41 (0) 44 8645111
Fax +41 (0) 44 8645005
info.vi@biotronik.com
www.biotronik.com

© 2018 BIOTRONIK AG – All rights reserved.
Specifications are subject to modification, revision and improvement.

 **BIOTRONIK**
excellence for life